

Adding Cytokine into Cells

A procedure for refreshing cell culture medium and adding cytokines, followed by microscopic inspection and returning the cells to the incubator.

Equipment

- CO2 Incubator (HERACELL VIOS 160i)
- TC Hood (Tissue culture hood for sterile work)
- Aspiration system (Vacuum line with glass pipette for medium removal)
- Pipette controller (Used with serological pipettes)
- Micropipette (Used for adding cytokines)
- Inverted microscope (EVOS system for cell imaging)

Reagents

- 70% ethanol (Used for hand disinfection)
- Fresh medium (Pre-warmed medium for cell culture)
- Cytokine(s)

Procedure

1. Disinfect hands thoroughly with 70% ethanol.
2. Carefully remove the required cell culture plat from the incubator, and trans port the plate into TC hood
3. Inside the hood, gently aspirate and discard the old medium from the culture plate without disturbing the cells.
4. Using a sterile pipette, add the appropriate volume of pre-warmed fresh mediu m to the plate.
5. Add cytokine(s) to the medium at the desired working concentration.
6. Examine the cells under an inverted microscope to assess cell health, density , and morphology.
7. Place the plate back into the incubator at the correct conditions.