

CRISPR Delivery

A procedure for delivering CRISPR components into 293T cells by mixing three reagents and adding the mixture to a cell culture dish.

Equipment

- 1.5 mL EP tube (sterile)
- pipette
- pipette tips
- 10 cm dish (containing 293T cells)

Reagents

- reagent 1
- reagent 2
- reagent 3
- 293T cells (In a 10 cm dish)

Procedure

1. Add reagent 1 into a sterile 1.5 mL EP tube.
2. Add reagent 2 into the EP tube.
3. Add reagent 3 into the EP tube.
4. Mix well.
5. Incubate at room temperature for 20 min.
6. Add the mixture dropwise into a 10 cm dish of 293T cells (~70-80% confluency)
7. Gently rock the dish forward and backward to mix.